

KERRACEL 3M GELLING FIBER DRESSING 10x10 cm



Product Code	59588
Unit of sale	box of 10 pcs
Minimum order	1
Type	Medical device
Class	II B

RDM (NSIS)	2569173
CND	M040404
PARAF	983199047
GTIN	04054596952683

Description

3M™ KERRACEL™

3M™ code: CWL1166

Description: Kerracel™ gelling fiber dressing 10x10 cm

Kerracel™ gelling fibre dressing is a versatile carboxymethylcellulose (CMC) primary dressing for moderately exuding wounds. Made of 100% carboxymethylcellulose (CMC) fibres, Kerracel™ dressing absorbs and retains high levels of exudate, forming a gel which conforms to the wound bed and supports a moist wound healing environment. As fluid is absorbed vertically and locked away, peri-wound skin is protected from maceration. The unstitched design aids microcontouring to the wound bed, minimising the space where bacteria can proliferate. Conformability aids ease of application, whilst its ability to retain integrity when wet helps to facilitate integral removal.