

3M COMPLY STERIGAGE STEAM CHEMICAL INTEGRATOR 1243A



Product Code	35959
Unit of sale	box of 1000 pcs
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
GTIN	50707387469073

Description

3M™ COMPLY™ STERIGAGE™ CHEMICAL INTEGRATOR - box of 1,000

3M™ Comply™ (SteriGage™) Chemical Integrator 1243A is a Class 5 Integrating Integrator that meet the performance requirements of ANSI/AAMI/ISO 11140-1:2005.

- moving front ink technology provides easy to read "Accept or Reject" result
- can be used in all 250-275°F steam sterilization cycles
- response correlates results correlate with a Biological

Indicator at three time/temperature relationships under ideal steam sterilization conditions.