

TEGADERM 3M HIGH INTEGRITY ALGINATE DRESSING 5x5 cm



Product Code 59591

Unit of sale box of 10 pcs

Minimum order 1

Type Medical device

Class II B

RDM (NSIS) 2056586

CND M04040201

GTIN 04054596862074

Description

3M™ TEGADERM™ ALGINATE DRESSINGS

3M™ code: 90110

Tegaderm™ high integrity alginate dressing 5x5 cm

3M™ Tegaderm™ High Integrity Alginate Dressing helps to maintain a moist wound healing environment in medium to heavily exuding wounds.

The dressing is highly conformable, soft, absorbent, sterile, primary wound dressing.

In reaction with wound exudate, it forms a gelatinous mass which maintains a moist healing environment.

The high integrity design helps to provide structure to the dressing upon removal.

The dressing is also intended to help control minor bleeding and support autolytic debridement.

It requires an appropriate secondary dressing.

Compatible with 3M cover dressings.

Indicated for:

- Pressure, venous and neuropathic (diabetic) ulcers
- Superficial wound, cuts and abrasions
- Donor sites
- Post operative wounds
- Trauma wounds
- Other dermal lesions