

TEGADERM 3M HIGH INTEGRITY ALGinate DRESSING TAPE 30.4x2 cm



Product Code	59593
Unit of sale	box of 5 pcs
Minimum order	1
Type	Medical device
Class	II B

RDM (NSIS)	2056608
CND	M04040201
GTIN	04054596862104

Description

3M™ TEGADERM™ ALGinate DRESSINGS

3M™ code: 90120
Tegaderm™ high integrity alginate dressing rope 30.4x2 cm

3M™ Tegaderm™ High Integrity Alginate Dressing helps to maintain a moist wound healing environment in medium to heavily exuding wounds.
The dressing is highly conformable, soft, absorbent, sterile, primary wound dressing.
In reaction with wound exudate, it forms a gelatinous mass which maintains a moist healing environment.
The high integrity design helps to provide structure to the dressing upon removal.
The dressing is also intended to help control minor bleeding and support autolytic debridement.
It requires an appropriate secondary dressing.
Compatible with 3M cover dressings.

- Indicated for:
- Pressure, venous and neuropathic (diabetic) ulcers
 - Superficial wound, cuts and abrasions
 - Donor sites
 - Post operative wounds
 - Trauma wounds
 - Other dermal lesions