

SCOTCHCAST 3M 5 cm x 3.65 m - green



Product Code	59611
Unit of sale	box of 10 pcs
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2002562
CND	M03050202
GTIN	50707387394818

Description

3M™ Scotchcast™ PLUS Casting Tape

- 3M Code: 82002
- Colour: Green
- Dimensions: 5 cm x 3.65 m

A lightweight, strong and durable casting tape that combines the benefits of a fiberglass casting tape with the handling ease of plaster.

- water-activated resin lets you control working time according to application needs
- high and early strength allows for weight bearing in 20 minutes and provides a durable cast with less chance of breakdown and the need for costly repairs Or replacements
- lightweight for easy mobility
- stretches to conform well to contours for a more comfortable, better fitting cast
- porous weave allows the skin to breathe