

SCOTCHCAST 3M 7.5 cm x 3.65 m - blue



|               |                |
|---------------|----------------|
| Product Code  | 59623          |
| Unit of sale  | box of 10 pcs  |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | I              |

|            |                |
|------------|----------------|
| RDM (NSIS) | 2002559        |
| CND        | M03050202      |
| GTIN       | 50707387394870 |

Description

3M™ Scotchcast™ PLUS Casting Tape

- 3M Code: 82003
- Colour: Blue
- Dimensions: 7.5 cm x 3.65 m

A lightweight, strong and durable casting tape that combines the benefits of a fiberglass casting tape with the handling ease of plaster.

- water-activated resin lets you control working time according to application needs
- high and early strength allows for weight bearing in 20 minutes and provides a durable cast with less chance of breakdown and the need for costly repairs Or replacements
- lightweight for easy mobility
- stretches to conform well to contours for a more comfortable, better fitting cast
- porous weave allows the skin to breathe