

URINE CONTAINER 120 ml - sterile



Product Code	25963
--------------	-------

Unit of sale	box of 250 pcs
--------------	-------------------

Minimum order	1
---------------	---

Type	IVD - Professional Use
------	---------------------------

Class	
-------	--

UK-REP	No
--------	----

CH-REP	No
--------	----

RDM (NSIS)	2518689
------------	---------

CND	W05010203
-----	-----------

EAN/UPC	8023279259636
---------	---------------

PARAF	976327092
-------	-----------

Description

Graduated urine containers made in polypropylene with red screw cap and frosted label.

Individually wrapped

Sterile

Capacity: 120 ml

Quantity: Box of 250 pcs.

One multilanguage instruction manual in each carton: GB, FR, IT, ES, PT.