

SCOTCHCAST 3M 10 cm x 3.65 m - blue



Product Code	59633
Unit of sale	box of 10 pcs
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2002560
CND	M03050202
GTIN	50707387394948

Description

3M™ Scotchcast™ PLUS Casting Tape

- 3M Code: 82004
- Colour: Blue
- Dimensions: 10 cm x 3.65 m

A lightweight, strong and durable casting tape that combines the benefits of a fiberglass casting tape with the handling ease of plaster.

- water-activated resin lets you control working time according to application needs
- high and early strength allows for weight bearing in 20 minutes and provides a durable cast with less chance of breakdown and the need for costly repairs Or replacements
- lightweight for easy mobility
- stretches to conform well to contours for a more comfortable, better fitting cast
- porous weave allows the skin to breathe