

SOFTCAST 3M 7.5 cm x 3.65 m - blue



Product Code	59673
Unit of sale	box of 10 pcs
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	1984635
CND	M03050202
GTIN	50707387463873

Description

3M™ Soft-Cast™ Casting Tape

- 3M Code: 82103
- Colour: Blue
- Dimensions: 7.5 cm x 3.65 m

A fiberglass casting tape that does not set to a completely rigid state, but remains flexible and soft even when it is completely cured.

- the soft, pliable cast edges improve patient comfort
- can be used in combination with rigid fiberglass splints, as soft edges provide additional comfort

Suggested Applications:

- sports Medicine: thumb, wrist and ankle sprains
- pediatric Orthopedics: serial casting for club foot treatment
- general Orthopedics: secondary casting, hybrid casting, corsets
- occupational Therapy: removable splints