

URINE VACUUM TUBE 10 ml - sterile



Product Code	25975
Unit of sale	box of 100 pcs
Minimum order	1
Type	IVD - Professional Use
Class	

UK-REP	No
CH-REP	No

RDM (NSIS)	1935692
CND	W050301020102
EAN/UPC	8023279259759

Description

Sterile urine vacuum tube - 10 ml

Material: Polypropylene

Individually wrapped urine test tube with paper label in single box.
Graduated 0.5, 1, 2.5, 5, 10 ml.
Multilingual box: GB, FR, IT, ES.

Vacuum tube yellow cap without preservative.