

# STOCKINET 3M 5 cm x 22.8 m



|               |                |
|---------------|----------------|
| Product Code  | 59701          |
| Unit of sale  | 1 pc           |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | I              |

|            |                |
|------------|----------------|
| RDM (NSIS) | 1949415        |
| CND        | M03030201      |
| GTIN       | 50707387788341 |

## Description

3M™ Synthetic Cast Stockinet™

- 3M Code: MS02
- STOCKINET™ 5 cm x 22.8 m

3M™ Stockinet is a synthetic tubular knit made of polyester fibers, used as protection for the patient's skin under synthetic plaster braces.  
It is not a sterile product but can be steam sterilized in an autoclave.

Water repellent, making it suitable for use with Soft Cast as a waterproof cast.  
Suitable for use with all casting materials, including plaster and synthetic.