

# UNDER-CAST PADDING 3M 15 cm x 2.7 m



Product Code	59718
Unit of sale	box of 12 pcs
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2233115
CND	M010201
GTIN	28711428073953

## Description

### 3M™ Synthetic Under-Cast Padding

- 3M Code: MW06
- Dimensions: 15 cm x 2.7 m

Porous and highly conformable non-woven polyester padding. Strong, yet easy to tear.

- synthetic cast padding
- non-woven polyester padding which is porous, and highly conformable
- easy to tear
- suitable for use with all casting materials