

NASOGASTRIC SECUREMENT DEVICE 3M 4.5x4.7 cm



| | |
|---------------|----------------|
| Product Code | 59725 |
| Unit of sale | box of 50 pcs |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|------------|----------------|
| RDM (NSIS) | 2391891 |
| PARAF | 964436341 |
| GTIN | 30707387792382 |

Description

3M™ NASOGASTRIC SECUREMENT DEVICE

- 3M™ code: 1500NG
- Nasogastric device 4.5x4.7 cm
- Box of 50

The 3M™ Nasogastric Securement Device is designed to provide reliable securement to help keep nasal feeding tubes and other tubes delivering therapy via the nose in place. It offers an easy-to-use solution that supports pressure injury prevention practices and minimizes the risk of medical adhesive-related skin injury (MARSİ).

Key Features:

- More Reliable Securement: Adhesive stays in place in challenging locations for up to 48 hours.
- Skin Performance: Designed to help reduce the risk of pressure injury by allowing the nasogastric tube to hang freely.
- Skin Performance: Designed to deliver the securement power needed while minimizing skin damage.
- Easy to Use: Easy to apply and remove; allows for standardization of application technique.
- Transparent: Provides visibility of the skin and nasal feeding tubes, aiding in assessment.
- Non-Adhesive Tab: The tab at the end of the tail is designed to assist in removal.