

NASOGASTRIC SECUREMENT DEVICE 3M 5.7x7,3 cm



Product Code	59726
Unit of sale	box of 50 pcs
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2391896
CND	M04010299
PARAF	964378780
GTIN	30707387792375

Description

3M™ NASOGASTRIC SECUREMENT DEVICE

- 3M™ code: 1500NG
- Nasogastric device 4.5x4.7 cm
- Box of 50

The 3M™ Nasogastric Securement Device is designed to provide reliable securement to help keep nasal feeding tubes and other tubes delivering therapy via the nose in place. It offers an easy-to-use solution that supports pressure injury prevention practices and minimizes the risk of medical adhesive-related skin injury (MARSi).

Key Features:

- More Reliable Securement: Adhesive stays in place in challenging locations for up to 48 hours.
- Skin Performance: Designed to help reduce the risk of pressure injury by allowing the nasogastric tube to hang freely.
- Skin Performance: Designed to deliver the securement power needed while minimizing skin damage.
- Easy to Use: Easy to apply and remove; allows for standardization of application technique.
- Transparent: Provides visibility of the skin and nasal feeding tubes, aiding in assessment.
- Non-Adhesive Tab: The tab at the end of the tail is designed to assist in removal.