

# 24 HOURS URINE BOTTLE 2500 ml



|               |                        |
|---------------|------------------------|
| Product Code  | 25980                  |
| Unit of sale  | box of 30 pcs          |
| Minimum order | 1                      |
| Type          | IVD - Professional Use |
| Class         | A                      |

|        |    |
|--------|----|
| UK-REP | No |
| CH-REP | No |

|            |               |
|------------|---------------|
| RDM (NSIS) | 2518682       |
| CND        | W05010203     |
| EAN/UPC    | 8023279259803 |
| PARAF      | 976327104     |

## Description

24 h graduated urine containers made in polyethylene.

Bottle 2500 ml

One multilanguage instruction in each carton: GB, FR, IT, ES, PT.