

2 LINE LABELLER MD20 10+10 DIGITS



Product Code	35985
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8034049916123
PARAF	990652822

Description

<div>2 LINE LABELLER MD20 10+10 DIGITS</div></div>

<div><div>Ergonomic and lightweight labeller made of high-quality durable plastic. With its advanced 10+10 digit format, the system instantly prints and applies labels containing essential data including: LOT number, sterilization date, expiry date, cycle number, operator, product reference and process indicator. The improved roll loading and faster setup for the operator reduce the risks of printing errors.</div></div>

<div><div>A COMPLETE TRACEABILITY OF THE STERILIZATION CYCLE</div><div>This medical labelling system is engineered to bring precision, safety, and efficiency to the sterile processing workflow. At a glance, it confirms that the instrument has been properly processed through sterilization: the labeller incorporates new print fonts that make it easier to read sterilization and expiration dates, all with the aim of preventing operator errors during data entry.</div></div>

<div><div>Made in Italy.</div></div>

</div>