

# URINE CONTAINER 120 ml in single box - sterile



Product Code	25987
Unit of sale	box of 100 pcs
Minimum order	1
Type	IVD - Professional Use
Class	

UK-REP	No
CH-REP	No

RDM (NSIS)	2518688
CND	W05010203
EAN/UPC	8023279259872
PARAF	976327080
GTIN	08023279259872

## Description

Graduated and sterile urine container with red cap - 120 ml.

Individually wrapped sterile containers with paper label and screw cap in single box.

Multilanguage box: GB, FR, IT, ES.

Material: Polypropylene