

URINE TEST TUBE 12 ml in single box - sterile



Product Code	25990
Unit of sale	box of 100 pcs
Minimum order	1
Type	IVD - Professional Use
Class	

UK-REP	No
CH-REP	No

RDM (NSIS)	2243403
CND	W050301020102
EAN/UPC	8023279259902
PARAF	976327256
GTIN	08023279259902

Description

Sterile urine test tube - 12 ml.

Individually wrapped urine test tube with paper label in single box.

Graduated 0.5, 1, 2.5, 5, 10 ml.

Material: Polypropylene

Multilanguage box: GB, FR, IT, ES.