

# 24 HOURS URINE TANK 2500 ml in single box



Product Code	25995
Unit of sale	box of 27 pcs
Minimum order	1
Type	IVD - Professional Use
Class	A

UK-REP	No
CH-REP	No

RDM (NSIS)	2518687
CND	W05010203
EAN/UPC	8023279259957
PARAF	976327270
GTIN	08023279259957

## Description

Graduated 24 h urine tank.

Polyethylene urine tank for 24 hour urine collection, in single box.

Multilanguage box: GB, FR, IT, ES.