

IMPRESSION TRAYS - set of 10 solid



| | |
|---------------|----------------|
| Product Code | 60045 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|--------|-----|
| UK-REP | Yes |
| CH-REP | Yes |

| | |
|------------|----------------|
| RDM (NSIS) | 2740290 |
| CND | V0402 |
| EAN/UPC | 8023279600452 |
| GTIN | 08023279600452 |
| GMDN | 42893 |

Description

Set of 10 stainless steel solid impression trays (5 lower + 5 upper)