

# FLAT REELS 200 m x 200 mm AUTOCLAVE



Product Code	36005
Unit of sale	box of 4 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes

RDM (NSIS)	2437114
CND	S01020101
GTIN	18033844700552
GMDN	13735

## Description

FLAT REELS FOR STERILIZATION BY STEAM OR ETO - 200 m

External size: 200 m x 200 mm

Sterilization reels produced with the best certified materials: high grammage medical paper (60 g/m<sup>2</sup>) and a double layer of polyethylene/polypropylene film (55 g/m<sup>2</sup>), in compliance with MDR (EU) 2017/745 specially designed to ensure high antibacterial protection for medical sterilization.

They guarantee a constant seal during sterilization; extreme ease of opening without fibre release or tearing; proper maintenance of sterility.

Medical Device Class I.

Indicators provide clear colour changes when exposed to sterilization process:

- steam: light blue to green
- ethylene oxide: pink to orange
- formaldehyde: green to yellow.