



Product Code	56015
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2799006
CND	V0899
EAN/UPC	6974580870237

Description

GIMA CAMP BED - foldable - green

Foldable and lightweight camp bed made of high strength aluminium alloy and 600D polyester. It is mainly used for army and battle field.

Supplied with carrying bag.

Multilingual manual: GB, FR, IT, ES, PT, DE, NL, DK, SE, NO, FI, LT, LV, EE, PL, CZ, SK, HU, SI, HR, RO, BG, GR, Arabic

Technical Specifications

Open size cm: 196x75x43
 Folded size cm: 95x20x11
 Weight: 6.4 kg
 Load: kg 159 kg