

CARTRIDGE SYRINGE - 2.2 cc



Product Code	60195
--------------	-------

Unit of sale	1 pc
--------------	------

Minimum order	1
---------------	---

Type	Medical device
------	----------------

Class	I
-------	---

UK-REP	No
--------	----

CH-REP	No
--------	----

RDM (NSIS)	2835657
------------	---------

CND	A020203
-----	---------

EAN/UPC	8023279601954
---------	---------------

GMDN	35969
------	-------

Description

Stainless steel Cartridge Syringe - 2.2 cc with thread in millimeter with aspiration.