

# SELF-SEAL POUCHES 60x110 mm



Product Code	36030
Unit of sale	box of 10000 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes

RDM (NSIS)	2180053
CND	S01010101
GTIN	18033844701047
GMDN	13735

## Description

### SELF-SEAL POUCHES FOR STERILIZATION BY STEAM OR ETO

External size: 60 x 110 mm

Minimum order: 10,000 pcs (Packed in box of 200 pouches)

Self-seal pouches provide very strong and equally divided sealings (multiline).  
The pouches have an adhesive band with the innovative system "finger-lift" which allows an easier use.  
Duo-check system: double sterilization indicator, internal and external.  
It ensures high control in the sterilization process.

Indicators provide clear colour changes when exposed to sterilization process:

- steam: light blue to green
- ethylene oxide: pink to orange
- formaldehyde: green to yellow.