

SELF-SEAL POUCHES 90x230 mm



Product Code	36032
Unit of sale	box of 3000 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes

RDM (NSIS)	2466025
CND	S01010101
GTIN	18033844701054
GMDN	13735

Description

SELF-SEAL POUCHES FOR STERILIZATION BY STEAM OR ETO

External size: 90 x 230 mm

Minimum order: 3,000 pcs (Packed in box of 200 pouches)

Self-seal pouches provide very strong and equally divided sealings (multiline). The pouches have an adhesive band with the innovative system "finger-lift" which allows an easier use. Duo-check system: double sterilization indicator, internal and external. It ensures high control in the sterilization process.

Indicators provide clear colour changes when exposed to sterilization process:

- steam: light blue to green
- ethylene oxide: pink to orange
- formaldehyde: green to yellow.