



Product Code	26049
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

EAN/UPC 7045430035362

## Description

### LAERDAL FEMALE MULTI - VEINOUS IV TRAINING KIT

Lifelike arm reproduction with multi-vein system designed for peripheral intravenous therapy.

Kit includes: arm, replacement skin and multi vein system, simulated blood, blood bag with tubing and connector, clamp and hook, 5 syringes, manikin lubricant, carrying case and directions for use.

Length 73.7 cm

Weight 1.79 kg