



Product Code	21605
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

EAN/UPC	8023279216059
---------	---------------

Description

JACKET WITH STUDS 60% COTTON - 40% POLYESTER

- Size: XL
- Colour: white
- UNISEX
- Five hidden press studs
- Short sleeves
- 2 front and 1 breast pockets

Professional and stylish line with all sizes for men and women, made of 60% high quality cotton and 40% polyester.
 Stock in 5 colours, other 22 colours available in 60 days.
 White jacket 200 g/m2, coloured jacket 185 g/m2.

Technical Specifications

Composition: 60% cotton - 40% polyester
 Weight: 200 g/m2 (white), 185 g/m2 (coloured)
 Yarn count: warp: Ne 20 ± 2 weft: Ne 18 ± 2
 Construction: warp: 40 ± 3 weft: 19 ± 3
 Breaking strength - warp: 100 kg ±3%
 Breaking strength - weft: 39 kg ±3% (white) 49 kg ±3% (coloured)
 Fastness to washing: 4
 Treatments and finishing: crease resist
 Dimensional stability - washing at 90°C: 3%