

NEEDLE-LESS VALVE - sterile



Product Code	26059
Unit of sale	box of 500 pcs
Minimum order	1
Type	Medical device
Class	I S

UK-REP	No
CH-REP	No

RDM (NSIS)	956617
CND	A07050202
EAN/UPC	8033395208906

Description

NEEDLE-LESS VALVE - sterile - box of 500

Device for I.V. infusion and/or withdrawal, needle-less closed circuit.
For protection of vascular access and of all infusion systems equipped with male luer-lock standard connections.
This closed and self-sealing protection system allows to infuse or inject fluid, make withdrawals eliminating the risk of accidental puncture to operators and reducing the risk of contamination to patients.

Latex free.
PVC free.