

Y INFUSION SET - polybag



Product Code	26062
Unit of sale	box of 500 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	2343580
CND	A03010101
EAN/UPC	8032764150358

Description

Y INFUSION SET - POLYBAG - box of 500 (20 bags of 25)

Infusion set includes: ABS spike with air vent, 6 cm PVC chamber, 150 cm PVC tube, ABS roller, ABS luer lock terminal without needle, latex free injection site.
Sterilized by Ethylene oxide.