

# WHITE COAT - cotton - unisex press studs size 48



Product Code	26084
Unit of sale	1 pc
Minimum order	1
Type	PPE
Class	I

UK-REP	No
CH-REP	No
EAN/UPC	8023279260847
PARAF	977811126
GTIN	08023279260847

## Description

Doctor's white coats 100% cotton

- Unisex, with press studs
- Size
  - Europe: 48
  - USA: M
  - France: 44

Professional and stylish line with all sizes for men and women, made of high quality 100% cotton. Doctor's white coats have 2 side pockets, 1 breast pocket, half belt and back slit. Woman and man models have buttons, unisex model has press studs and elasticated cuffs. Oeko-Tex® 100 certified. Weight 185 g/m2.