

WHITE COAT - cotton - unisex press studs size 50



Product Code	26085
Unit of sale	1 pc
Minimum order	1
Type	PPE
Class	I

UK-REP	No
CH-REP	No
EAN/UPC	8023279260854
PARAF	977811138
GTIN	08023279260854

Description

Doctor's white coats 100% cotton

- Unisex, with press studs
- Size
 - Europe: 50
 - USA: L
 - France: 46

Professional and stylish line with all sizes for men and women, made of high quality 100% cotton. Doctor's white coats have 2 side pockets, 1 breast pocket, half belt and back slit. Woman and man models have buttons, unisex model has press studs and elasticated cuffs. Oeko-Tex® 100 certified. Weight 185 g/m2.