



Product Code	26233
Unit of sale	1 pair
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes
RDM (NSIS)	2837004
CND	V9099
EAN/UPC	8023279262339
PARAF	909847612
GTIN	08023279262339
GMDN	13576

Description

GIMA PROFESSIONAL CLOGS MADE IN ITALY

- UPPER WITHOUT PORES
- WITH STRAP

- active arch support to stimulate circulation
- antistatic material directly injected in the polymeric blend
- perforated sole for breathability and transpiration
- antislip grooved sole safe on wet floor
- machine washable, sterilizable at 134°C
- available also with strap for extra security
- individually packed in carton box
- latex free