

# "ARGENTA" TORCH - metal - silver



Product Code	25633
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
--------	----

RDM (NSIS)	2755589
CND	V9099
EAN/UPC	8023279256338
PARAF	903120424
GMDN	12347

## Description

Professional penlights, also suitable for diagnostic, with bulb or white LED with quick clip shifter mechanism or push button (25634). Argenta and Elegance made of anodized aluminium, Nova in topquality ABS plastic housing. Supplied with 2x1.5V AAA batteries.

Argenta torch  
Light: bulb  
Light intensity: 4.5 lux  
Material: metal  
Supplied with batteries: 2 x AAA  
Colour: silver  
ON/OFF by: Pocket clip

Made in Germany.