



Product Code	26332
Unit of sale	1 pair
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2836983
CND	V9099
EAN/UPC	8023279263329
PARAF	903598732
GTIN	08023279263329
GMDN	13576

## Description

GIMA PROFESSIONAL CLOGS MADE IN ITALY

• UPPER WITHOUT PORES

- active arch support to stimulate circulation
- antistatic material directly injected in the polymeric blend
- perforated sole for breathability and transpiration
- antislip grooved sole safe on wet floor
- machine washable, sterilizable at 134°C
- available also with strap for extra security
- individually packed in carton box
- latex free