



| | |
|---------------|----------------|
| Product Code | 25637 |
| Unit of sale | box of 50 pcs |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|------------|---------------|
| RDM (NSIS) | 2219438 |
| CND | T020604 |
| EAN/UPC | 8030413002959 |

Description

MASK-OP® with wrapping visor

- lacets
- light blue

Liquid splash resistant and anti-fogging surgical mask with antistatic visor.

- 4 layers: three made of non-allergenic woven nonwoven and one of polyethylene net
- nuisance-free and adaptable nose-clip
- filtering efficiency (B.F.E.) =99%, according to EN 14683:2019+AC:2019
- resistance to breath <30,0 Pa/cm², according to EN 14683:2019+AC:2019
- latex-free
- non sterile
- liquid splash resistance: resistant at pressure of 16 kPa, according to ISO 22609, related to EN 14683:2019+AC:2019
- anti-fogging internal stripe made of expanded polyethylene
- Bioburden = 30 CFU/g, according to EN 14683:2019+AC:2019

Medical Device Class I Regulation EU 2017/745.

Multilingual box: GB, FR, IT, ES, PT, DE, GR, PL, RU, Arabic.

Made in Italy.