



Product Code	26373
Unit of sale	1 pair
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2836996
CND	V9099
EAN/UPC	8023279263732
PARAF	903599037
GTIN	08023279263732
GMDN	13576

Description

GIMA PROFESSIONAL CLOGS MADE IN ITALY

- UPPER WITHOUT PORES

- active arch support to stimulate circulation
- antistatic material directly injected in the polymeric blend
- perforated sole for breathability and transpiration
- antislip grooved sole safe on wet floor
- machine washable, sterilizable at 134°C
- available also with strap for extra security
- individually packed in carton box
- latex free