



Product Code	26474
Unit of sale	1 pair
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes
RDM (NSIS)	2837016
CND	V9099
EAN/UPC	8023279264746
GTIN	08023279264746
GMDN	13576

Description

GIMA PROFESSIONAL CLOGS MADE IN ITALY

- UPPER WITH PORES
- WITH STRAP

- active arch support to stimulate circulation
- antistatic material directly injected in the polymeric blend
- perforated sole for breathability and transpiration
- antislip grooved sole safe on wet floor
- machine washable, sterilizable at 134°C
- available also with strap for extra security
- individually packed in carton box
- latex free

Technical Specifications

- En ISO 20347 A-E