

# GOLD SCISSORS STRAIGHT BLUNT/SHARP - 14.5 cm



Product Code	26500
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes

RDM (NSIS)	128220
CND	L010499
EAN/UPC	8023279265002
PARAF	950833842
GTIN	08023279265002
GMDN	38727

## Description

### GOLD LINE - TUNGSTEN CARBIDE INSTRUMENTS

A top quality theatre use Tungsten Carbide Instruments made of stainless steel AISI 410 (AISI 420 for scissors) with very hard tips made in Germany to help surgeons perform precision surgeries with best results.

Highly reliable for hardness and corrosion resistance.

Instructions: GB, FR, IT, ES, PT, DE, GR, Arabic.