



Product Code	26591
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes

RDM (NSIS)	2837815
CND	Y091299
EAN/UPC	8023279265910
PARAF	903162410
GTIN	08023279265910
GMDN	40504

Description

Stainless Steel AISI 304 holloware with minimum thickness 0.6 mm
Male vertical urinal
Size: Ø 114 x 220 mm - 1.5 l

Technical Specifications

- Autoclavable 121°C
- AISI 304