

S/S INSTRUMENT TRAY - 380X304X50 mm



Product Code	26597
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2740258
CND	V0402
EAN/UPC	8023279265972
GTIN	08023279265972
GMDN	12143

Description

Stainless Steel AISI 304 holloware with minimum thickness 0.6 mm
Instrument tray - angled edges - without cover
Size: 380 x 304 x h 50 mm

Technical Specifications

- Autoclavable 121°C
- AISI 304
- Without cover