

S/S INSTRUMENT TRAY - 223X126X45 mm



Product Code	26599
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2740258
CND	V0402
EAN/UPC	8023279265996
PARAF	903162459
GTIN	08023279265996
GMDN	12143

Description

Stainless Steel AISI 304 holloware with minimum thickness 0.6 mm
Instrument tray - round edges - without cover
Size: 223 x 126 x h 45 mm

- Autoclavable 121°C
- AISI 304
- Without cover