

# STOCKINETTE LEGGINGS 35x120 cm - large



Product Code	23660
Unit of sale	box of 40 pcs
Minimum order	1
Type	Medical device
Class	I S

UK-REP	No
CH-REP	No

RDM (NSIS)	2120004
CND	T0299
EAN/UPC	8023279236606
GMDN	35531

## Description

STOCKINETTE LEGGING 35x120 cm - large - sterile - box of 40

For orthopedics and gynaecology

Multilingual instructions and label: GB, FR, IT, ES, PT, DE, NL, DK, SE, NO, FI, LT, LV, CE, PL, CZ, SK, HU, SI, HR, RO, BG, GR, Arabic

Made in Italy.