

S/S INSTRUMENT TRAY - 306X196X50 mm



Product Code	26603
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2740258
CND	V0402
EAN/UPC	8023279266030
PARAF	903162485
GTIN	08023279266030
GMDN	12143

Description

Stainless Steel AISI 304 holloware with minimum thickness 0.6 mm
Instrument tray - angled edges - without cover
Size: 306 x 196 x h 50 mm

- Autoclavable 121°C
- AISI 304
- Without cover