

S/S DRESSING JAR 0.5 l with lid - diam.106x66 mm



Product Code	26623
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2740284
CND	V0402
EAN/UPC	8023279266238
PARAF	903162802
GTIN	08023279266238
GMDN	42893

Description

Stainless Steel AISI 304 holloware with minimum thickness 0.6 mm
Dressing jar - with lid
Size: Ø 102 x 65 mm - 0.5 l

- Autoclavable 121°C
- AISI 304
- With cover