

S/S DRESSING JAR 2 l with lid - diam.127x162 mm



Product Code	26625
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2740284
CND	V0402
EAN/UPC	8023279266252
PARAF	903162840
GTIN	08023279266252
GMDN	42893

Description

Stainless Steel AISI 304 holloware with minimum thickness 0.6 mm
Dressing jar - with lid
Size: Ø 127 x 161 mm - 2 l

- Autoclavable 121°C
- AISI 304
- With cover