

# S/S FORCEPS JAR - diam. 55x180 mm



|               |                |
|---------------|----------------|
| Product Code  | 26627          |
| Unit of sale  | 1 pc           |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | I              |

|        |     |
|--------|-----|
| UK-REP | Yes |
| CH-REP | Yes |

|            |                |
|------------|----------------|
| RDM (NSIS) | 2740308        |
| CND        | V0402          |
| EAN/UPC    | 8023279266276  |
| PARAF      | 903162865      |
| GTIN       | 08023279266276 |
| GMDN       | 42893          |

## Description

Stainless Steel AISI 304 holloware with minimum thickness 0.6 mm

Forceps jar  
Size: Ø 55 x h 180 mm  
Autoclavable 121°C