

# S/S THERMOMETER JAR - diam. 33x80 mm



Product Code	26628
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes

RDM (NSIS)	2831158
CND	V0399
EAN/UPC	8023279266283
PARAF	903176651
GTIN	08023279266283
GMDN	42893

## Description

Stainless Steel AISI 304 holloware with minimum thickness 0.6 mm

Thermometer jar  
Size: Ø 33 x h 80 mm  
Autoclavable 121°C