

GALLIPOT/LOTION BOWL 100 mm - plastic - graduated 300 ml



Product Code	26632
Unit of sale	1 pc
Minimum order	10
Type	Medical device
Class	I

UK-REP	Not required
CH-REP	Not required

RDM (NSIS)	2656170
CND	V0402
EAN/UPC	8023279266320
PARAF	903162939

Description

PLASTIC HOLLOWARE - Autoclavable hundreds of times

Gallipot/lotion bowl
Size: 100 x h 45 mm
Graduation: 300 ml

For the protection and delivery of sterile instruments to operating rooms. For the administration of procedures, holding cleansing solutions and the administration of medicines and fluids.

4 good reasons to choose plastic line:

- 50% lighter than stainless steel
- autoclavable: 1000 cycles guaranteed
- attractive modern design
- cheaper than equivalent stainless steel

Made of polypropylene material, easy to clean, resistant to most healthcare chemicals. Machine, ultrasonic or manual cleaning are permissible. Autoclave sterilization at 134.6°C (+2°C) for a minimum of 2 minutes. EO Ethylene Oxide or Hydrogen Peroxide sterilization both at max. 65°C.

Made in the UK.