

INSTRUMENT TRAY 300X250X52 mm - plastic



| | |
|---------------|----------------|
| Product Code | 26639 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|--------|-----|
| UK-REP | Yes |
| CH-REP | No |

| | |
|------------|----------------|
| RDM (NSIS) | 2627977 |
| CND | V0402 |
| EAN/UPC | 8023279266399 |
| PARAF | 903162992 |
| GTIN | 08023279266399 |
| GMDN | 12143 |

Description

PLASTIC HOLLOWWARE - Autoclavable hundreds of times

Instrument Tray - Solid Ribbed Base

Size mm: 300x250x52

Graduation: -

For the protection and delivery of sterile instruments to operating rooms. For the administration of procedures, holding cleansing solutions and the administration of medicines and fluids.

4 good reasons to choose plastic line:

- 50% lighter than stainless steel
- autoclavable: 1000 cycles guaranteed
- attractive modern design
- cheaper than equivalent stainless steel

Made of polypropylene material, easy to clean, resistant to most healthcare chemicals. Machine, ultrasonic or manual cleaning are permissible. Autoclave sterilization at 134.6°C (+2°C) for a minimum of 2 minutes. EO Ethylene Oxide or Hydrogen Peroxide sterilization both at max. 65°C.

Made in the UK.