



Product Code	26642
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2739353
CND	V0402
EAN/UPC	8023279266429
GTIN	08023279266429
GMDN	42893

## Description

Stainless Steel holloware, autoclavable 121°C, with 0.5 mm thickness.

Economic kidney dish - with curved edges - tray 8"

- Size: 207 x 98 x 39 mm
- Capacity: 400 ml

## Technical Specifications

- Autoclavable 121°C
- AISI 304